
medisana®

Non-contact & Multi-functional Infrared Thermometer

TM 320

- GB Unit and controls**
- 1 Display
 - 2 On/Off-button
 - 3 Sensor cover
 - 4 Measurement result
 - 5 Forehead measurement
 - 6 Measurement-button
 - 7 Battery compartment
 - 8 Sound symbol
 - 9 Measurement unit (°C or °F)
 - 10 Ear measurement
 - 11 Memory-button
 - 12 Sensor
 - 13 Battery symbol
 - 14 Object measurement
 - 15 Memory display

Explanation of symbols

This instruction manual belongs to this device. It contains important information about starting up and operation. Read the instruction manual thoroughly. Non-observance of these instructions can result in serious injury or damage to the device.

WARNING
These warning notes must be observed to prevent any injury to the user.

CAUTION
These notes must be observed to prevent any damage to the device.

NOTE
These notes give you useful additional information on the installation or operation.

Classification: Type BF applied part

Temperature range

Keep dry

IP22 Protection classification against solid objects and water

LOT LOT number

Manufacturer

Date of manufacture

48624_TM 320_manual_GB_202205_L11

GB Infrared Multi Functional Thermometer

Manual - Please read carefully!

Safety Information
IMPORTANT INFORMATION! RETAIN FOR FUTURE USE!
Read the instruction manual carefully before using this device, especially the safety instructions, and keep the instruction manual for future use. Should you give this device to another person, it is vital that you also pass on these instructions for use.

INTENDED USE:
Thank you for purchasing our Non-contact & Multi-functional Infrared Thermometer TM 320. This Infrared Thermometer is a non-sterile, reusable, handheld device. It can be used by consumers in home care environment and doctors in clinic as reference. It is intended for measuring human body temperature of three months and above of people by detecting infrared heat from the center of the forehead and in the ear.

- Only use the unit for the specific purpose described in the instruction manual. Any misuse will void the warranty.
- The unit is not designed for commercial use.
- Using this thermometer is no substitute for consulting a doctor.
- The thermometer and the patient must remain in a stable environment for at least 30 minutes before measuring the temperature.
- This device is not designed to be used by persons (including children) with limited physical, sensory or mental abilities, or by persons with insufficient experience and/or knowledge, unless under observation by a person responsible for their safety, or unless they have been instructed in the use of the device.
- Children must be supervised to ensure that they do not play with the device.
- The temperature in the forehead and temple area differs from the internal temperature, which is taken orally or rectally. If the measurement therefore does not match the patient's own perception or is unusually low, repeat the measurement every 15 minutes. As a reference, you can also measure the internal body temperature using a conventional oral or rectal thermometer.
- A control measurement using a conventional thermometer is recommended in the following cases:
 1. If the reading is surprisingly low
 2. For new-born infants up to 100 days old
 3. For children under three years of age who have a weakened immune system or who react unusually in the presence or absence of fever
 4. If you are using the thermometer for the first time or want to familiarise yourself with it
- If a fault occurs or the unit is damaged, do not attempt to repair the unit yourself, as this will void the warranty. Contact your dealer and only have repairs carried out by authorized service partners.
- Do not modify and / or open the device. It contains small parts which might be swallowed by children.
- Do not use the thermometer if it is damaged.
- Do not use the thermometer in damp environments. Do not immerse the thermometer in water or other liquids. It is not watertight.
- This device comprises sensitive components and must be treated with caution. Observe the storage and operating conditions described in the «Technical Specifications» section.
- Do not use this device close to strong electromagnetic fields such as mobile telephones or radio installations.
- Keep a minimum distance of 3.3 m from such devices when using this device.
- Do not shake or tap the thermometer. Do not drop it.
- Protect the device from direct sunlight, extremely high or low temperatures, dirt and dust.
- If the thermometer is not going to be used for a long period, remove the battery.

WARNING - The measurement results obtained by this thermometer do not represent a diagnosis and do not replace the consultation of a doctor, who should always be consulted in case of doubt.

Items supplied and packaging
Please check first of all that the device is complete and is not damaged in any way. In case of doubt, do not use it and contact your dealer or your service centre. The following parts are included:
 • 1 medisana Non-contact & Multi-functional Infrared Thermometer TM 320
 • 2 Batteries (AAA type, LR03) 1.5V
 • 1 Instruction manual

The packaging can be reused or recycled. Please dispose properly of any packaging material no longer required. If you notice any transport damage during unpacking, please contact your dealer without delay.

WARNING - Please ensure that polythene packing is kept away from the reach of children!
Risk of suffocation!

Taking a temperature reading - quickly and accurately!
The medisana Non-contact & Multi-functional Infrared Thermometer allows precise, contact-free measurement of body temperature (at a distance of approximately 0-5 mm). The temperature is taken above the eyebrow between the forehead and the temple.
The value of the energy radiated there is measured and converted to a temperature value. The reading is displayed after 2 seconds. This enables you, for example, to take babies' temperature without disturbing their sleep.

Additional function: Apart from measuring body temperature, the Infrared multifunctional thermometer can also be used to measure the temperature of objects (e.g. fluids, surfaces, ambient temperature). The temperature range for this function is 0°C- 93.2°C (32°F-199.7°F).

Operation - Getting started
Remove/ insert battery: If the battery symbol 15 appears in the LC display after a certain period of use, you must replace the battery by a new 1.5V battery (AAA type, LR03). Removal: Press the cover of battery compartment 7, take it off and take out the battery. Insertion: Put in a new battery. Pay attention to correct polarity (please refer to instruction inside the battery compartment 7). Put the battery compartment 7 cover back on.

BATTERY SAFETY INFORMATION

- Do not disassemble batteries!
- Clean the battery and device contacts if necessary before putting in the batteries!
- Remove discharged batteries from the device immediately!
- Increased risk of leakage, avoid contact with skin, eyes and mucous membranes! If battery acid comes in contact with any of these parts, rinse the affected area with copious amounts of fresh water and seek medical attention immediately!
- If a battery has been swallowed seek medical attention immediately!
- Insert the battery correctly, observing the polarity!
- Keep the battery compartment well sealed!
- Remove the battery from the device if it is not going to be used for an extended period!
- Keep batteries out of childrens reach!
- Do not attempt to recharge batteries! There is a danger of explosion!
- Do not short circuit! There is a danger of explosion!
- Do not throw into a fire! There is a danger of explosion!
- Keep unused batteries in their packaging away from metal objects in order to prevent short circuiting!
- Do not throw used batteries into the household refuse; put them in a hazardous waste container or take them to a battery collection point, at the shop where they were purchased!

NOTES
The device will automatically switch off if there is no input for approx. 30 seconds.

Before use
Acclimatise the thermometer to the temperature of the room: For an accurate reading of a patient's temperature, their body must first acclimatise to the temperature of the room where the measurement will be made. The same applies to the Multi-functional Non-contact Infrared Thermometer TM 320. When taken somewhere else, it automatically adjusts itself to the new environment. The thermometer takes a certain amount of time to equalise its own temperature to that of the room. This depends on the difference between the two temperatures and can take more than 30 minutes. Precise readings can only be taken once this process is completed.

Measuring forehead temperature


NOTES For the most accurate readings:

- Do not take a baby's temperature during or directly after breastfeeding.
- Patients should not drink or eat directly before or during a measurement, and they should avoid any physical exertion.
- Clean the area where you are taking the temperature and remove any dirt, hairs or sweat before moving the sensor 4 to the measuring position.
- Do not remove the thermometer from the measuring position until it has beeped to signal the end of the reading.
- If the sensor becomes dirty during or after a reading, clean it with a cotton bud dipped in cleaning alcohol before putting the thermometer away.
- Always take the temperature in the same place, because the readings will vary in different environments.


The thermometer's safety and accuracy has been inspected and confirmed, provided it is used in accordance with these instructions.

1. Put on sensor cover 3 and then press the On/Off-button 2 to switch on the thermometer. Now the last measured temperature value appears along with the corresponding symbol (for forehead temperature and for in-ear temperature). After approx. 3 seconds the device is ready for measurement - this will be indicated by a short beep.
2. To begin the measurement, point the sensor 4 at the middle of the patient's forehead at a distance of 0-5 mm. If the patient's forehead has hair, sweat or dirt on it, remove this first so that the reading is as accurate as possible.
3. Press the Measurement-button 6. You will see the LCD is flashing in the display 1.
4. After approx. 1 second a long beep can be heard to indicate that the measurement is finished and the measured temperature will be shown.
5. Remove the sensor from the measuring position. If the result is <34°C, it shows LO with red backlight. When the result is 34°C-37.1°C, it shows green backlight. When the result is 37.2°C-38.1°C, it shows orange backlight. When the result is 38.2°C-43.0°C, it shows red backlight. When the result is >43.0°C, it shows and HI with red backlight.

After approximate 5 seconds the device is ready for another measurement. The device will automatically switch off if there is no input for approx. 30 seconds. You may also switch the device off manually by pressing the On/Off-button 2.

- Measuring the temperature of in-ear**
1. To measure the temperature in your ear, first remove the sensor cover 3.
 2. Then press the On/Off-button 2 to switch on the thermometer. Now the last measured temperature value appears along with the corresponding symbol (for forehead temperature and for in-ear temperature). After approx. 3 seconds the device is ready for measurement - this will be indicated by a short beep.
 3. Stretch your ear canal by pulling your ear backwards and upwards and carefully insert the sensor into the ear. For children under the age of 1 year, the sensor should be inserted into the ear canal, while the child is lying flat (e.g. on a bed) with the head slightly turned, so that the ear points upward.
 4. Press the Measurement-button 6. You will see the LCD is flashing in the display 1.
 5. After approx. 1 second a long beep can be heard to indicate that the measurement is finished and the measured temperature will be shown.
 6. If the result is <34°C, it shows LO with red backlight. When the result is 34°C-37.1°C, it shows green backlight. When the result is 37.2°C-38.1°C, it shows orange backlight. When the result is 38.2°C-43.0°C, it shows red backlight. When the result is >43.0°C, it shows and HI with red backlight.

- Measuring the temperature of an object**
1. Put on sensor cover 3 and then press the On/Off-button 2 for at least 6 seconds. Now the last measured temperature value appears along with the corresponding symbol (for forehead temperature and for in-ear temperature). After approx. 3 seconds the device is ready for measurement - this will be indicated by a short beep.
 2. To begin the measurement, point the sensor 4 at the targeted object.
 3. Press the MEM-button 11. You will see the LCD is flashing in the display 1.
 4. After approx. 1 second a long beep can be heard to indicate that the measurement is finished and the measured temperature will be shown.
 5. Remove the sensor from the measuring position. The measurement result will show with white backlight. If the result is <0°C, it shows LO. When the result is >93.2°C, it shows HI.

After approximate 5 seconds the device is ready for another measurement. The device will automatically switch off if there is no input for approx. 30 seconds and will automatically switched to body measurement mode. You may also switch the device off manually by pressing the On/Off-button 2.

Change the measurement unit
You can measure the temperature in °C or in °F. Please proceed as follows to select the unit of measurement:
 1. The device has to be switched off.
 2. Long press the MEM-button 11 until you see the symbol «°C» or «°F» flashing.
 3. Press the MEM-button 11 again to change the unit of measurement.
 4. Press the On/Off-button 2 to save the selected measurement unit. Or wait for 8 seconds to save the selected measurement unit and the device is ready for the next measurement.

Recalling stored data
The thermometer can store up to 9 results. You can call up your stored data with a short press of the MEM-button 11, without having previously switched on the thermometer with the On/Off-button 2. The symbol appears in the display, together with "On".
If measurement values are stored in memory, first the memory slot number «1» appears and afterwards the stored temperature with the corresponding symbol (for body temperature and object temperature). By further pressing the MEM-button 11, you can call up more stored values (memory slot «2» etc.). Totally there are 9 memory slots available, which may be recalled by pressing the MEM-button 11 repeatedly.

Audio settings
You may activate or deactivate the acoustic signals of the device. To adjust, press the On/Off-button 2 to switch on the thermometer. The last measured temperature value appears along with the corresponding symbol (for forehead temperature and for in-ear temperature).
Now you may adjust the signal settings by pressing the MEM-button 11 until the signals are activated (symbol:) or deactivated (symbol will be disappeared). Afterwards, press the On/Off-button 2 to save the setting and switch off the device.

Error messages

Screen display	Reasons	Solutions
HI	When the target object temperature is higher than measurement range which the ear temperature is more than 43.0°C, the LCD screen will show "HI" indication.	1. When the measurement probe is not properly placed in the ear canal or if the measurement distance is too far during the measurement process, the measurement result may be low. 2. When measurement probe is dirty, the measurement value may be low; we should use a cotton swab to dip alcohol to wash the measurement probe properly.
Lo	When the target object temperature is lower than measurement range which the ear temperature is less than 34°C, the LCD screen will show "Lo" indication.	
Er.H	This product has an upper limit of operating temperature of 40°C. When the ambient temperature exceeds this temperature point, the LCD screen will display an error message "Er.H".	When you operate this product, the environment temperature cannot be greater than 40°C.
Er.L	This product has a lower limit of operating temperature of 10°C. When the ambient temperature exceeds this temperature point, the LCD screen will display an error message "Er.L".	When you operate this product, the environment temperature can not be lower than 10°C.

When the operation is wrong, if the sound is turned on, then the sound will read the value and remind at this moment. (A short beep will be heard.)

Troubleshooting

Phenomenons	Reasons	Solutions
When the power is turned on, the screen can not be displayed.	The battery is exhausted.	Replace the new battery.
	The battery polarity is wrong.	The battery polarity is the same as the battery case.
The measurement temperature is low.	The measurement position is not correct.	Measure the temperature correctly according to the instructions.
	There are dirt dust stuff in the sensor or ear canal.	Please clear the dirt before the measurement.
Big temperature fluctuations with continuous measurement.	The measurement interval is too small.	The interval for each measurement should be above 10 seconds.

Cleaning and maintenance

- The sensor 4 is the most delicate part of the device. Protect it from dirt and damage. Clean the sensor with a cotton bud dipped in warm water or 70% isopropyl alcohol.
- Do not immerse the thermometer in water or any other liquid. It is not watertight.
- To clean the thermometer, use a soft cloth slightly dampened with mild soapy water or 70% isopropyl alcohol. Do not let water get into the unit. Only use the unit again when completely dry.
- Never use aggressive cleaning agents, thinners, benzine or tough brushes.
- Take care not to scratch the lens of the sensor, the display or the housing.
- Do not use or store the unit in an environment where the ambient temperature or humidity is too high or too low (see operating conditions in the "Technical specifications" section), in direct sunshine, in the vicinity of static electricity or in a dusty environment.
- Remove the batteries if you do not wish to use the unit for an extended period of time. Otherwise there is a danger that the batteries may leak.
- If the thermometer is used as described in these instructions, no regular maintenance or calibration is required.

Disposal
This product must not be disposed of together with domestic waste. All users are obliged to hand in all electrical or electronic devices, regardless of whether or not they contain toxic substances, at a municipal or commercial collection point so that they can be disposed of in an environmentally acceptable manner.
Please remove the battery before disposing of the device/unit. Do not dispose of old batteries with your household waste, but at a battery collection station at a recycling site or in a shop. Consult your municipal authority or your supplier for information about disposal.

Technical specifications

Name and model:	medisana Non-contact & Multi-functional Infrared Thermometer TM 320
Display system:	digital display (LCD)
Power supply:	2 x batteries (AAA type, LR03) 1.5V
Expected service life:	5 years
Measuring range:	Body mode 34-43°C (93.2- 109.4°F); Object mode 0°C-93.2°C (32°F-199.7°F)
Accuracy of measurement (laboratory readings):	± 0.2°C, 35.0°C-42.0°C (± 0.4°F, 95.0°F-107.9°F), other temperature ± 0.3°C
Resolution:	0.1°C (0.1°F)
Memory capacity:	up to 9 readings
Operating temperature:	10.0°C-40.0°C (50.0°F-104.0°F); rel. humidity 15-90% RH
Storage/transport temperature:	-25°C to +55°C (-13°F-131°F), rel. humidity 0-93%
Dimensions:	approx. 31 x 175 x 72 mm
Weight:	approx. 77 g
Article No.:	48624
EAN Code:	40155 88 48624 1

The CE-mark of this device refers to the EC-directive 93/42/EEC MDD. The unit is designed for home use only. For commercial use, a metrological inspection must be carried out once every year. This calibration will be charged for and can be carried out by an appropriate authority or an authorised service centre - in accordance with the regulations for users of medical equipment. Electromagnetic compatibility: The device complies with the EN 60601-1-2 standard for electromagnetic compatibility. Details on this measurement data can be found at the end of this manual.

In accordance with our policy of continual product improvement, we reserve the right to make technical and optical changes without notice.

Warranty and repair terms
Please contact your dealer or the service centre in case of a claim under the warranty. If you have to return the unit, please enclose a copy of your receipt and state what the defect is. The following warranty terms apply:
 1. The warranty period for medisana products is three years from date of purchase. In case of a warranty claim, the date of purchase has to be proven by means of the sales receipt or invoice.
 2. Defects in material or workmanship will be removed free of charge within the warranty period.
 3. Repairs under warranty do not extend the warranty period either for the unit or for the replacement parts.
 4. The following is excluded under the warranty:
 a. All damage which has arisen due to improper treatment, e.g. nonobservance of the user instructions.
 b. All damage which is due to repairs or tampering by the customer or unauthorised third parties.
 c. Damage which has arisen during transport from the manufacturer to the consumer or during transport to the service centre.
 d. Accessories which are subject to normal wear and tear.
 5. Liability for direct or indirect consequential losses caused by the unit are excluded even if the damage to the unit is accepted as a warranty claim.

Shenzhen Pango Medical Electronics Co., Ltd.
2-4 Floor, No.5 Shanzhuang Rd., Xikeng Village, Henggang Town, Longgang District, Shenzhen, Guangdong, China.

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